

Introduction: Fully automated Microcarrier based Suspension Cultivation

The increased demand for cells that more closely mimic in vivo functions require novel cell culture approaches. We developed a complete automated 3D cell culture system that we used to streamline the cell culture process of human Neural Stem Cells (hNSC) (Fig. 1):

GEMs (Global Eucaryotic Microcarriers) are microcarriers uniquely designed for cell culture. The BioLevigator supports the complete cell culture process from inoculation to harvesting. The 3D CellHOST automates parallel scale-up of different cultures at a medium-scale level.




Fig. 1: The 3D CellHOST System. The novel fully automated 3D scale-up system (3D CellHOST™) employs magnetic microcarriers (GEM™), an integrated incubator (BioLevigator™) and a liquid handling Workstation (STARlet™) for automated sampling and media changes. The technology is based on levitation of microcarriers in rotating tubes (LeviTubes™) and provides optimal access to nutrients and oxygen under low shear forces. The integrated BioLevigator ensures homogenisation of the cells / GEM suspension and provides all functions of a cell culture incubator.

Results: Fully automated 3D CellHOST mediated Expansion of Human Neural Stem Cells

1. Evaluation of 3D CellHOST mediated Expansion of human ES cell-derived Neural Stem Cultures (hESNSC)

✓ Cultivated cells: hESNSC



Fig. 2: Long-term self renewing human ES cell derived Neural Stem Cells (It-hESNSC). Recently we derived a neural stem cell population from H9.2 human Embryonic Stem Cells (hESC), which exhibit extensive self-renewal, clonogenicity, and stable neuro-genesis. Proliferating cells express high levels of the neural stem cell markers Nestin, Sox2, Sox1 and Pax6. Under optimal conditions only few cells differentiate spontaneously (Koch et al., 2009).

✓ 3D CellHOST cultivated hESNSC yielded higher cell numbers per area than their 2D counterparts




Fig. 3: Growth analysis of It-hESNSC after 3D CellHOST-mediated expansion. 2×10^4 / cm² It-hESNSCs were inoculated and subsequently cultured for 9 days either with 0.5 ml (3D Standard) or 2 ml (3D Scale-Up) of GEM slurry per culture tube. (A) Cultures were automatically sampled at day 1, 3, 6, 9 of cultivation and total cell number per tube were calculated from those samples. (B) Inoculation of 24×10^6 It-hESNSC (equivalent to two confluent 10 cm culture dishes) with 2 ml of GEM slurry per LeviTube and subsequent cultivation for 9 days within the 3D CellHOST (4 LeviTubes) yielded a 15-fold increase in cell number (equivalent to thirty 10 cm dishes). (C) The cell yield per cm² in the 3D CellHOST was compared to conventional 2D maintenance cultures. The determined cell number per cm² after 3D CellHOST mediated expansion was comparable higher as indicated.

2. Characterization of It-hESNSCs after fully automated 3D CellHOST-mediated expansion

✓ It-hESNSC maintain their typical morphology and metabolic activity. Their proliferation rate under scale-up conditions is enhanced.




Fig. 4: It-hESNSC maintain their typical rosette-like growth-pattern and metabolic activity. (A, B) It-hESNSC cultured for 9 days on laminin coated GEMs maintained their typical rosette-like morphology (C) like in a conventional 2D culture centred around tight-junction protein ZO-1⁺ lumens. (D) Glucose consumption and lactate production rates of It-hESNSCs cultured by the 3D CellHOST system are comparable to those of manually cultured cells.




Fig. 5: BrdU incorporation analysis in It-hESNSC during 3D CellHOST-mediated expansion. It-hESNSCs were cultured under 3D standard, 3D scale-up and 2D conditions as indicated. Cell cultures were analysed at day 3 and day 9 as displayed. It-hESNSCs were pulsed for 2 h with 10 μM BrdU. The analysis revealed comparable percentages at day two of cultivation for all cultures. Determination of BrdU positive cells at day 9 revealed more BrdU positive cells under scale-up conditions (E: 38.1 %) compared to 3D standard (D: 24.3 %) and 2D control (F: 20.5 %).

✓ It-hESNSC remain undifferentiated after fully automated 3D CellHOST-mediated expansion.




Fig. 6: Expression of the pan-neural markers Sox2 and Nestin. (A-F) According to immunocytochemical analysis It-hESNSC grown on laminin coated GEMs for 9 days maintain neural stem cell marker expression as indicated. Virtually all cells were positive for the pan-neural markers Nestin and Sox2 (C, F) like their 2D cultured counterparts.




Fig. 7: Abundance of differentiated cells. Immunocytochemical analysis of It-hESNSC, grown on laminin coated GEMs for 9 days showed no expression of (A-C) neuronal (TUJ1) or (D-F) glial (GFAP) markers after 3D CellHOST-mediated expansion.

Conclusion and Outlook

Conclusion: Fully automated 3D CellHOST-mediated cultivation of It-hESNSC:

- Facilitated higher cell yields per cm² (Figure 3).
- Was shown to maintain typical morphological characteristics and metabolic activity of It-hESNSC, thus demonstrating its biological compatibility (Figure 4).
- Enhanced the proliferation abilities of high-density cultures (Figure 5).
- Supported undifferentiated growth of It-hESNSC (Figure 6 and Figure 7).

Outlook : Automated 3D cell culture opens the path towards using cells as reagents.

- Uncoupled cell culture from downstream assays.
- From freezer to assay in one step.
- Automated 3D CellHOST mediated culture could facilitate parallel scale-up of pluripotent Stem Cells and their differentiated progeny.



Fig. 8: Flowchart of the cells as reagent model. Automated 3-D cell culture eliminates the need for seed-split cycles and streamlines the cell culture process. It also enables the uncoupling of cell culture from downstream assays and supports the direct use of cells from the freezer into an assay in a single step.